

SHIELD^{me}™

DISPOSABLE ISOLATION GOWN

SDM-DGWN

CE



SHIELD^{me}™ NON-SURGICAL ISOLATION GOWN (WITH BELT) is intended to protect health care personnel and patients from the transfer of microorganisms, body fluids (Fluid repellent) and particulate material. Single-use and disposable. Latex Free.

Gown made of non-woven fabrics using spun bonded polypropylene spunbond non-woven, 40 gsm, 100 % new chip, coated with polyethylene film (28 gsm PP + 12 gsm PE), and latex free material with back closure, waste and neckties, featured sleeves closed with a belt on the right side. Liquid resistant due to polypropylene spun bond non-woven coated with Polyethylene film

CE-label according to EU Directive for Medical equipment 93/42/EEC, Annex VII, class 1, non-sterile.

SPECIFICATIONS:

- Color** : Blue/Yellow
- Gown Length** : 110 Cm Regular and 115 Cm Ex-Large
- Storage** : Store under ambient temperature and dry condition, keep away from Direct Sunlight Durability unopened 5 years
- Waste Management** : Polybag, combustion in incineration plants only creates carbon dioxide and water, which directly forms part of the natural circulation. Cardboard is combusted or recycled
- Packing** : 10 pcs/Bag : 5 Bags/Carton (50 Gown / Carton)

AL ARIQ EQUIPMENT LLC

Phone: 24502925 E-mail: info@alariqoman.com

Muscat Sahar Salalah

99885172 99889374 99330886

www.alariqoman.com