

# SHIELD<sup>me</sup>™

CE

## DISPOSABLE PROTECTIVE SUIT

SME-DCL01



**SHIELD<sup>me</sup>™ DISPOSABLE PROTECTIVE SUIT** is a disposable anti-microbial protective clothing suitable for production, outdoor protection, sentry boxes, laboratories, dust-free workshops and peripheral personnel of hospitals and infectious contaminated area.



- With isolated and protective effects, it is also suitable for clinical personnel who come into contact with the blood, body fluids and secretions of potentially infectious patients when working in a non-sterile environment.
- Independent packaging design, it is easy to use and fits comfortably.
- Made of thermoplastic polyurethane elastomer TUP antibacterial film material, which is safe and reliable.
- Disposable hooded isolation suit has elastic wrists, ankles and waist.

### SPECIFICATION AND PACKING

**Sterilized Protective Clothing**

**Fabric weight:** 50g

With Hood without Boot

Water Proof, Oil-proof, Acid-Proof Liquid

Splash-proof

**Available Sizes:** M 160-165cm, L 170cm,  
XL 175cm, XXL 180cm XXXL 185cm

**Packing:** 1 pcs/Nylon Bag/50 pcs/Carton  
50x30x40cm

**Weight:** 10kgs/Carton

**AL ARIQ EQUIPMENT LLC**

Phone: 24502925 E-mail: info@alariqoman.com

Muscat

Sohar

Salalah

99885172

99889374

99330886

www.alariqoman.com